METHODOLOGY Open Access

A 3D human co-culture to model neuron-astrocyte interactions in tauopathies



Kevin L. Batenburg¹, Claudia Sestito^{2,3}, Paulien Cornelissen-Steijger^{3,4}, Jan R. T. van Weering^{1,5}, Leo S. Price², Vivi M. Heine^{3,4} and Wiep Scheper^{1,5*}

Abstract

Background Intraneuronal tau aggregation is the major pathological hallmark of neurodegenerative tauopathies. It is now generally acknowledged that tau aggregation also affects astrocytes in a cell non-autonomous manner. However, mechanisms involved are unclear, partly because of the lack of models that reflect the situation in the human tauopathy brain. To accurately model neuron-astrocyte interaction in tauopathies, there is a need for a model that contains both human neurons and human astrocytes, intraneuronal tau pathology and mimics the three-dimensional architecture of the brain.

Results Here we established a novel 100–200 µm thick 3D human neuron/astrocyte co-culture model of tau pathology, comprising homogenous populations of hiPSC-derived neurons and primary human astrocytes in microwell format. Using confocal, electron and live microscopy, we validate the procedures by showing that neurons in the 3D co-culture form pre- and postsynapses and display spontaneous calcium transients within 4 weeks. Astrocytes in the 3D co-culture display bipolar and stellate morphologies with extensive processes that ensheath neuronal somas, spatially align with axons and dendrites and can be found perisynaptically. The complex morphology of astrocytes and the interaction with neurons in the 3D co-culture mirrors that in the human brain, indicating the model's potential to study physiological and pathological neuron-astrocyte interaction in vitro. Finally, we successfully implemented a methodology to introduce seed-independent intraneuronal tau aggregation in the 3D co-culture, enabling study of neuron-astrocyte interaction in early tau pathogenesis.

Conclusions Altogether, these data provide proof-of-concept for the utility of this rapid, miniaturized, and standardized 3D model for cell type-specific manipulations, such as the intraneuronal pathology that is associated with neurodegenerative disorders.

Keywords hiPSC-derived neurons, Astrocytes, Neuron-astrocyte interaction, 3D co-culture, Tau aggregation

Wiep Scheper

w.scheper@amsterdamumc.nl

¹ Department of Functional Genomics, Center for Neurogenomics and Cognitive Research, Vrije Universiteit Amsterdam, Amsterdam Neuroscience - Neurodegeneration, De Boelelaan 1085, 1081 HV Amsterdam, The Netherlands

² Crown Bioscience Netherlands B.V. (Formerly OcellO B.V.), Leiden, The Netherlands

⁴ Department of Complex Trait Genetics, Center for Neurogenomics and Cognitive Research, Vrije Universiteit Amsterdam, Amsterdam Neuroscience, De Boelelaan 1085, 1081 HV Amsterdam, The Netherlands ⁵ Amsterdam UMC location Vrije Universiteit Amsterdam, Department of Human Genetics, Amsterdam Neuroscience - Neurodegeneration, De Boelelaan 1085, 1081 HV Amsterdam, The Netherlands



© The Author(s) 2023. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/ficenses/by/4.0/. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

^{*}Correspondence:

³ Amsterdam UMC location Vrije Universiteit Amsterdam, Department of Child and Adolescent Psychiatry, Amsterdam Neuroscience, De Boelelaan 1085, 1081 HV Amsterdam, The Netherlands

Background

Tauopathies are a class of neurodegenerative diseases that are characterized by the progressive aggregation of tau protein and include Alzheimer's Disease (AD) and frontotemporal dementias (FTD) such as Pick's Disease, Corticobasal Degeneration and Progressive Supranuclear Palsy. Since tau is abundantly expressed in neurons where it predominantly accumulates in disease conditions, it disrupts neuronal physiology and leads to neurodegeneration in a cell-autonomous manner. In addition, astrocytes are affected by neuronal tau pathology and actively contribute to neurodegeneration by the loss of homeostatic function and gain of toxic function [1-3]. Rodent models have generated important insight into how cell (non-) autonomous disease pathways converge and eventually lead to neurodegeneration [3–5]. However, mouse models are relatively expensive and often only recapitulate limited aspects of human disease [6]. More importantly, many therapeutic targets derived from rodent studies have thus far failed in clinical trials, emphasizing the need to establish human models of tau aggregation.

Human induced pluripotent stem cell (hiPSC) technology has created the possibility to study disease mechanisms in human neurons and astrocytes. The physical and functional interaction between these cell types, for example at tripartite synapses, is critical for brain physiology and plays an important, yet not fully understood role in neurodegeneration [7]. Recently our group developed a two-dimensional (2D) hiPSC-derived neuron/astrocyte co-culture model for tau pathology, to study cell nonautonomous effects of intraneuronal tau pathology [8]. Although hiPSC-derived neural cultures contributed to our understanding of neuron and astrocyte physiology, they lack the natural interactions in the brain microenvironment. This is a limitation for astrocytes to acquire their in vivo-like morphology and thus the complexity by which they interact with neurons in a three-dimensional (3D) network [9, 10]. This led to the development of 3D neuron/astrocyte models that display increased cellular and spatial complexity [11, 12], but these systems typically depend on the self-organization into organoids that are complex, time-consuming, costly and may suffer from heterogeneity and limited throughput. Organoids are valuable assets to investigate neurodevelopmental disorders [13], but are thus far less well suited for the study of aging-associated neurodegenerative disease. Moreover, the robust induction of intraneuronal tau pathology has not been reported in 3D organoids.

Therefore, the aim of this study is to establish a physiologically relevant and robust model representing 3D neuron-astrocyte interactions. To this end, we present a step-by-step procedure to generate $100-200~\mu m$ -thick 3D microwell co-cultures of hiPSC-derived neurons and

human primary astrocytes. This system is miniaturized, compatible with cell type-specific manipulation as well as standard imaging techniques, and can be efficiently generated at reduced cost. We show that neurons in the 3D co-culture form a morphologically mature synapse pattern and display spontaneous calcium activity within 4 weeks. Moreover, like in the human brain, astrocytes display a complex morphology with extensive processes that enwrap neuronal somas, extend upon axons and dendrites, and localize perisynaptically. As a proof-of-concept for the cell type-specific introduction of pathology in neurons, we successfully established seed-independent intraneuronal tau pathology in the 3D co-culture to enable the study of the role of neuron-astrocyte interactions early in the pathogenesis of tauopathies. Altogether, we established a rapid and robust miniaturized model that mimics the 3D context of the human brain to study the interaction between neurons and astrocytes and how this is altered in tauopathies.

Results

Development of a 3D human neuron/astrocyte co-culture

We developed a protocol for the rapid generation of miniaturized and homogeneous 3D neuron/astrocyte co-cultures (Supplementary Document 1). Briefly, we employed a clonal hiPSC-line with doxycyclin-inducible Neurogenin 2 (Ngn2) [14], that was previously shown to facilitate the efficient generation of glutamatergic neurons from hiPSCs [15, 16]. Ngn2-hiPSCs were pre-differentiated to neural progenitors for two days prior to 3D co-culturing with human primary astrocytes. This timewindow of separate neuron/astrocyte culturing enables cell type-specific manipulations such as viral transduction. Next, neural progenitors were detached and mixed with human primary astrocytes in medium mixed with Geltrex extracellular matrix (ECM) and supplemented with NT3 and BDNF for neuronal maturation for 4 weeks (Fig. 1A). Since we observed proliferation of astrocytes in the 3D culture, the neuron/astrocyte ratio was optimized to prevent overgrowth and subsequent clumping of astrocytes. To this end, 30.000 neural progenitors and 5.000 human astrocytes were mixed per well: a cell ratio and density that resulted in a homogeneous distribution in the 3D co-culture, thus facilitating subsequent microscopy analyses. Also, we observed that a cell suspension with 50% Geltrex (v/v), but not lower concentrations, polymerized into sufficiently thick 3D co-cultures that remained stable for at least 4 weeks. Finally, miniaturization of the protocol enabled the maintenance of 50 µL 3D co-cultures that are shaped as a thin cylinder on the bottom of 96-well microplates. In this format the co-cultures are easy to handle and can be maintained by media refreshment once every week. In addition, since




Fig. 1 Development of a 3D human neuron/astrocyte co-culture. A Schematic representation of the protocol to generate the 3D human neuron/astrocyte co-culture. Ngn2-hiPSCs are plated on day -2 for neural induction in stem cell (TeSR) medium containing ROCK inhibitor (RI), doxycyclin (dox), NT3 and BDNF. One day after medium (DMEM) change on day -1, Ngn2-neural progenitors are detached on day 0 and mixed with human primary astrocytes in medium and Geltrex for neuronal maturation in medium (NB + B27) containing NT3 and BDNF with weekly media refreshments for 4 weeks. B 100 μm-thick maximum intensity confocal image projection of neurons and astrocytes in 3D co-culture. Immunostaining was performed for neurons (NeuroChrom), astrocytes (GFAP) and presynapses (SYP1). Single channels are shown as a heatmap corresponding to the focal plane of the fluorescence (bottom, red; top, blue). C 160 μm-thick maximum intensity confocal image projection of neurons and astrocytes in 3D co-culture. Immunostaining was performed for neurons (β-3-tubulin, green) and astrocytes (GFAP, magenta). Single channels are shown in greyscale and the merged image includes nuclei (DAPI, blue). A zoom of the boxed area indicated in the merged image is shown

the microwell plates with a high-clarity foil bottom are compatible with imaging, this enabled the immediate processing of the 3D co-cultures within the plate. Altogether, these characteristics of our 3D neuron/astrocyte platform are assets to increase the throughput and reduce handling and costs.

To investigate neuronal differentiation in 3D co-culture with human astrocytes, we performed immunofluorescent labeling using a combination of neuronal and astrocytic markers after 4 weeks of 3D co-culture. These data show that the 3D ECM permits the diffusion of antibodies at similar concentrations used for 2D cultures and

confirm the compatibility of the 3D system with confocal microscopy. Confocal imaging within the 100 μm -thick 3D matrix showed the presence of the pan-neuronal marker NeuroChrom and presynaptic synaptophysin 1 (SYP1) puncta (Fig. 1B) as well as abundant expression of β -3-tubulin (Fig. 1C) throughout the 3D co-culture, confirming the generation of neuronal extensions and presynapses. Importantly, GFAP-positive astrocytes with multiple long processes were observed dispersed throughout the 3D neuronal network (Fig. 1B, C), indicating the integration of astrocytes in the 3D co-cultures. These data demonstrate the utility of this protocol for the generation of 3D co-cultures of human neurons and astrocytes.

Morphologically complex astrocytes interact with neurons in 3D co-culture

Astrocytes in vivo display a complex morphology and contain numerous processes [17]. GFAP immunoreactivity may not label all astrocytes and their thinner processes [18], therefore we performed additional immunostaining of the glial membrane marker CD44 to investigate astrocyte morphology in more detail. In agreement with a previous study in post-mortem human brain [19] these data confirmed that not all CD44-positive astrocytes express GFAP (Fig. 2A). Also, these data demonstrate that extensive thin CD44-positive astrocytic processes are present in the 3D co-culture. Finally, and in accordance with previous studies [20, 21], these data show that human astrocytes in our 3D co-culture exhibited a bipolar or stellate morphology (Fig. 2A), akin to the morphology of fibrous astrocytes in situ [22] that is typically not found in 2D co-culture [11]. These data demonstrate that human astrocytes in the 3D co-culture obtain the morphological complexity similar to those in the human brain.

In the human brain, astrocytes contact neuronal cell bodies and synapses via their extensive processes [17]. Although immunostaining for CD44 labeled more astrocytes and thinner processes compared to GFAP (Fig. 2A), it also increased non-specific staining of larger spherical structures, likely derived from membranous debris. Therefore, to minimize non-specific staining for the assessment of neuron-astrocyte contact in the 3D co-cultures we performed double immunofluorescent staining of GFAP and NeuroChrom to label astrocytes and neurons, respectively. Confocal microscopy demonstrated that astrocytes are closely associated with neurons and often enwrap neuronal somas and neurites (Fig. 2B and Supplementary Video 1). In addition, the astrocytes were regularly found aligned with neuronal processes (Fig. 2C and Supplementary Video 1). To further assess astrocyte contact with axons and dendrites, we performed immunofluorescence staining and confocal microscopy for axons (SMI312) and dendrites (MAP2). First, this showed the differential localization of axonal and dendritic markers in neuronal extensions (Fig. 2D), demonstrating that hiPSC-derived neurons in the 3D co-culture acquire polarity of axons and dendrites within 4 weeks. In addition, this demonstrated that astrocytic GFAP-positive processes align in close contact with both axon and dendrite tracts. (Fig. 2D). Altogether, these data demonstrate that human astrocytes in the 3D co-culture display an in vivo-like morphology and interact with hiPSC-derived neurons.

Neurons form synapses and display spontaneous calcium transients in 3D co-culture

Synapses permit the specific activity and function of the neuronal network and the interaction between astrocytes and synapses is pivotal for synapse function [23]. We demonstrated the formation of presynaptic structures throughout the 3D co-culture (Fig. 1B). To further investigate the formation and contact between pre- and postsynapses in the 3D co-culture, double immunofluorescent staining and subsequent confocal microscopy using presynaptic SYP1 and postsynaptic density protein 95 (PSD95) was performed. SYP1 and PSD95 positive puncta colocalized onto axons (NF200; Fig. 3A), demonstrating that neurons in this 3D co-culture form a morphologically mature synapse pattern [24, 25]. To confirm the contact between pre- and postsynaptic compartments in the 3D co-culture at higher resolution we next optimized the processing of 3D co-cultures for transmission electron microscopy (EM) imaging. EM imaging showed presynaptic terminals containing synaptic vesicles connected to an aligned postsynaptic plasma membrane with an electron-dense postsynaptic density (Fig. 3B). This confirms the formation of pre- and postsynapses in the 3D co-culture and demonstrates the compatibility of this system with ultrastructural analysis. Interestingly, EM also demonstrated cellular processes containing intracellular glycogen granules, a characteristic of astrocytes [26], in contact with neuronal compartments (Fig. 3C). This included close contact with presynaptic terminals containing synaptic vesicles connected to an aligned postsynaptic plasma membrane, suggestive of a pre- and postsynaptic terminal in a tripartite synapse. These data demonstrate synapse formation and further confirm the interaction between astrocytes and neurons in the 3D human co-culture.

We next examined whether hiPSC-derived neurons in the 3D co-culture display calcium transients, indicative of neuronal activity [27]. Using live-cell imaging of Fluo-5 AM, a calcium-sensitive dye, spontaneous calcium transients were observed and recorded for 10 min (Fig. 3D and Supplementary Video 2). After live-cell imaging the




Fig. 2 Morphologically complex astrocytes interact with neurons. HiPSC-derived neurons and human primary astrocytes were co-cultured in 3D for 4 weeks, and astrocyte morphology and interaction with neurons was assessed by confocal microscopy. A 160 μm-thick maximum intensity confocal image projection of neurons and astrocytes in 3D co-culture. Immunostaining was performed for astrocytes (GFAP, magenta; CD44, green). Single channels are shown in greyscale. A zoom of the boxed area indicated in the merged image is shown. B 112 μm-thick maximum intensity confocal image projection of neurons and astrocytes in 3D co-culture. Immunostaining was performed for neurons (NeuroChrom, green) and astrocytes (GFAP, magenta). Single channels are shown in greyscale. A zoom of the boxed area indicated in the merged image is shown. See Supplementary Video 1 for a confocal z-stack. C 105 μm-thick maximum intensity confocal image projection of neurons and astrocytes in 3D co-culture. Immunostaining was performed for neurons (NeuroChrom, green) and astrocytes (GFAP, magenta). A merged image is shown, in which three numbered boxed areas are shown as zoomed areas. In the zoomed images, the arrows indicate astrocytes extending upon neuronal processes. D 78 μm-thick maximum intensity confocal image projection of axons, dendrites and astrocytes in a 3D co-culture. Immunostaining was performed for axons (SMI312, green), dendrites (MAP2, blue) and astrocytes (GFAP, magenta). Single channels are maximum intensity projections and shown in greyscale. The boxed area indicated in the MAP2 channel is shown as a merged image and the arrows indicate an astrocyte interacting with axons and dendrites

sample was fixed and processed for immunofluorescence to determine the cellular source of the 15 calcium traces (Supplementary Fig. 1A-C) using β -3-tubulin and GFAP to distinguish neurons and astrocytes, respectively (Fig. 3E). This showed that neurons are the primary source of the recorded calcium traces (12/15 cells, 2 cells

unidentified; Supplementary Fig. 1C). One astrocyte displayed a very slow rise in Fluo-5 AM intensity throughout the 10 min-record, suggesting the gradual internalization of the dye rather than astrocytic calcium activity. For a representative active neuron, the burst-to baseline calcium transients lasted 3–5 s (Fig. 3F), in accordance with




Fig. 3 Synapse formation and neuronal activity. Synapse formation (A-C) and calcium activity (D-F) was assessed in 4-week old 3D human neuron/astrocyte co-cultures A 12 µm-thick maximum intensity confocal image projection of axons, pre- and postsynapses in a 3D co-culture. Immunostaining was performed for axons (NF200, blue), pre- (SYP1, green) and postsynapses (PSD95, magenta). Single channels are shown in greyscale and zooms of the numbered boxed areas, indicated in the NF200 panel, are shown as merged. Intensity profiles of SYP1 (green) and PSD95 (magenta) were determined along the line segment indicated in merged zoom 3 to show co-localization of pre- and postsynaptic compartments. B Transmission electron micrograph of a pre- and postsynapse (pseudocolour in green and red, respectively). The boxed area indicated in the pseudocoloured image is shown as a greyscale image. C Transmission electron micrograph of an astrocyte (blue) in contact with a pre- and postsynaptic compartment (pseudocolour in green and red, respectively). Note that the astrocyte contains electron-dense glycogen granules. The boxed area indicated in the pseudocoloured image is shown as a greyscale image. A 4-week old 3D human neuron/astrocyte co-culture **D**, **F** was loaded with Fluo5-AM for calcium imaging and **E** subsequently labeled with β-3-tubulin and GFAP to assess the neuronal/ $astrocytic source of the calcium signal. \textbf{D} \ A \ 10-min \ maximum \ fluorescence intensity \ projection \ of \ Fluo 5-AM. The \ heatmap \ corresponds to the \ angle \ an$ calcium concentration. See Supplementary Video 2 for calcium signal over time. E After calcium imaging, immunofluorescent staining was performed for neurons (β-3-tubulin, green) and astrocytes (GFAP, magenta), and nuclei were visualized using DAPI (greyscale). F Calcium bursting in one cell indicated by the boxed area in D at three different time points (t=16, 20 and 24 s) and the corresponding immunofluorescent image of the boxed area in E. The graph shows the Fluo5-AM fluorescence intensity of this cell (cell 2 in Supplementary Fig. 1A-C) over 1 min, with the lowest and highest value set to 0 and 1 respectively, and the blue box corresponds to the timeframe between t = 16 and 24 s for which the bursting activity is shown. Please note that no calcium signal was observed in the astrocyte. See Supplementary Fig. 1 for the neuronal/astrocytic annotation of calcium traces

the timescales previously observed for human neurons cultured in 3D microfluidic platforms [20, 28]. These data demonstrate that the human neurons in the 3D co-culture form synapses and display neuronal activity.

Intraneuronal tau aggregation in 3D co-culture

In tauopathies, tau aggregates are present in neurons and contribute to neurodegeneration in a cell-autonomous and non-autonomous manner. To investigate

whether the 3D co-culture model can be employed for the study of cell type-specific pathology, we introduced tau pathology specifically in neurons. The introduction of tau containing two FTD-associated mutations P301L and S320F (FTDtau¹⁺²) previously enabled seed-independent tau aggregation in 3D rodent brain slice cultures [29]. Viral transduction of neural progenitors with FTDtau¹ (P301L) or FTDtau¹⁺² (P301L+S320F) was performed one day before 3D co-culturing with human astrocytes. Confocal microscopy showed that expression of FTDtau¹⁺² but not FTDtau¹ resulted in the formation of MC1-positive, pathological tau in MAP2positive neurons in the 3D co-culture within 4 weeks (Fig. 4A). Moreover, MC1-positive FTDtau¹⁺² was positive for the thiophene dye HS169 [30] (Fig. 4B), demonstrating the formation of β -sheet-rich tau aggregates.

The protocol also allows the introduction of differentially manipulated neuronal populations to investigate neuron-neuron interactions. As a proof-of-concept a mixture of neurons expressing FTDtau¹- and FTDtau¹+² with distinct fluorescent tags (EGFP and mCherry) was introduced in 3D co-culture with human astrocytes. These data show that FTDtau¹- and FTDtau¹+²-expressing neurons could be clearly distinguished and importantly only the latter population developed MC1-positive tau pathology (Fig. 4C), confirming that the phenotype of different tau variants (Fig. 4A) is reproduced in a mixed co-culture.

Tau aggregation in the 3D co-culture is independent on exogenous seeds and therefore models the effects of intraneuronal tau pathology. This makes the model potentially suitable for the study of cell non-autonomous disease mechanisms. Interestingly, interactions of GFAP-positive astrocytes with neurites and cell bodies of β -3-tubulin-positive neurons containing FTDtau $^{1+2}$ are present (Fig. 4D). This indicates that the interaction between neurons and astrocytes is intact. Altogether these data provide proof-of-concept for the utility of

the model to study early, cell non-autonomous disease mechanisms in the pathogenesis of tauopathies.

Discussion

Here we established and validated a novel 3D human neuron/astrocyte co-culture model of tau pathology. First, we report a step-by-step procedure to generate $100-200~\mu m$ -thick microwell co-cultures of hiPSC-derived neurons and primary human astrocytes. Using live, confocal as well as electron microscopy techniques, we demonstrate that the neurons in the 3D co-culture form pre- and postsynapses, display spontaneous calcium transients, and are in close contact with astrocytes. Second, by the introduction of two pathogenic tau mutations we are the first to demonstrate seed-independent tau aggregation in a human 3D co-culture.

We present a step-by-step protocol to generate miniaturized 3D co-cultures. These 3D co-cultures are more robust than traditional 3D organoids that may suffer from heterogeneity in size and cellular composition, limited throughput, and culturing and processing procedures that may be time-consuming or costly [31]. The 3D co-cultures presented here consist of homogeneous populations of hiPSC-derived neurons and human primary astrocytes and are compatible with standard imaging techniques. First, hiPSC-derived neural progenitors are mixed in a defined ratio with primary human astrocytes and subsequently embedded within an ECM, resulting in a homogeneous distribution and consistent density of cells. Second, since the ECM polymerizes at the bottom of 96-well microplates, it shapes as a thin cylinder that permits easy media refreshment and is structurally supported by the well sides. The 3D ECM remained stable for at least 4 weeks but preliminary data indicate that co-cultures are stable at 8 weeks, suggesting that they can also be used for studies that require prolonged culturing. Third, by utilizing 96-well plates we miniaturized the system and the high-clarity foil bottom ensured its compatibility with standard imaging techniques with high signal

(See figure on next page.)

Fig. 4 Seed-independent intraneuronal tau aggregation. Neural progenitors were transduced with mCherry-tagged FTDtau¹ or EGFP-tagged FTDtau¹⁺², mixed separately (**A, B, D**) or combined (**C**) with human primary astrocytes and tau pathology was assessed at week 4. **A** 12 μm-thick maximum intensity confocal image projections of 3D-cocultures expressing EGFP-tagged FTDtau¹ or FTDtau¹⁺² (green), immunostained for neurons (MAP2, blue) and a pathological conformation of tau (MC1, magenta). Single channels are shown in greyscale. A zoom of the boxed area indicated in the MC1 channel is shown as merged image and shows overlap of FTDtau¹⁺² and MC1 within a MAP2-positive neuron. **B** 11 μm-thick maximum intensity confocal image projection of 3D co-cultures expressing EGFP-tagged FTDtau¹ or FTDtau¹⁺², immunostained for a pathological conformation of tau (MC1), co-labelled for stacked beta-sheets (HS169). Nuclei were visualized using DAPI (blue) and are merged with FTDtau (grey). The HS169 and MC1 single channels are shown in greyscale and the merged image shows that only FTDtau¹⁺² (blue) overlaps with MC1 (magenta) and HS169 (green). **C** 10 μm-thick maximum intensity confocal image of a 3D neuron/astrocyte co-culture with neurons expressing either mCherry-tagged FTDtau¹ (blue) or EGFP-tagged FTDtau¹⁺² (green), immunostained for a pathological conformation of tau (MC1, magenta). A zoom of the boxed area indicated in the MC1 channel is shown as merged image and shows that only FTDtau¹⁺² overlaps with MC1. **D** 50 μm-thick maximum intensity confocal image projections of a 3D-coculture expressing EGFP-tagged FTDtau¹⁺² (green), immunostained for neurons (β-3-tubulin, blue) and astrocytes (GFAP, magenta). Single channels are shown in greyscale. A zoom of the boxed area indicated in the GFAP channel is shown as merged image. The arrows in the orthogonal views show an astrocytic extension in contact with a FTDtau¹⁺²-containing neuron




Fig. 4 (See legend on previous page.)

to noise. This format also permits the processing as well as (live) imaging of the 3D co-cultures directly in the plate, as opposed to organoids that typically require collection, embedding and sectioning. Specific 3D culture plates are commercially available, in which cells are either pipetted into microfluidic channels [20] or plated on top of a precast ECM matrix [32]. In contrast, in the current study standard cell culture plates are used, thus reducing technical and financial limitations to broadly implement the protocol.

Both the neurons and astrocytes present mature phenotypes within 4 weeks. Using immunofluorescence, live- and EM imaging techniques, we demonstrate that hiPSC-derived neurons in the 3D co-culture express typical neuronal markers, acquire polarity of axons and dendrites, form pre- and postsynapses and display calcium activity within 4 weeks. In addition, we show that astrocytes in the 3D co-culture are integrated in the neuronal network and display bipolar and stellate morphologies with extensive processes that are not typically found in 2D monoculture. This is consistent with previous studies that showed that rodent [33] as well as human [34] astrocytes adopt a more in vivo-like morphology when maintained in a 3D compared to 2D environment [10]. Also, rodent studies have previously demonstrated that morphological complexity of astrocytes is increased by co-culture with neurons [35, 36], largely mediated via connexins [37] and neuroligins [35] in a contact-dependent manner. Similarly, and in agreement with previous observations in an hiPSC-derived organoid system [11], we demonstrate that astrocytes interact with neurons in the 3D co-culture and extend their processes along axon and dendrite tracts. Importantly, astrocytes were found in close contact with neuronal cell bodies as well as with synapses, suggesting the formation of tripartite synapses. Indeed, such contacts are pivotal for astrocytes to control neuronal homeostasis in the brain, which includes the regulation of ion levels, storage and supply of energy, turnover of neurotransmitters as well the formation, function and elimination of synapses [23, 26, 38]. Altogether, the 3D co-culture offers a novel tool to study both cell types and the mechanisms by which they interact in a physiologically relevant 3D environment.

To our knowledge, we are the first to demonstrate intraneuronal tau pathology in a 3D human co-culture model in a seed-independent manner. The study of aging-related diseases such as sporadic tauopathies in hiPSC-derived neurons is challenging as they rarely match the transcriptional and functional maturity of human neurons in vivo [15, 39], and lack postnatal 4R tau isoforms even after 100 days in vitro [40, 41]. As is commonly observed in developing neurons, hiPSC-derived neurons often present with increased levels of tau phosphorylation, but do

not develop insoluble aggregates of endogenous tau [40, 42-45]. Altogether, these factors necessitate the overexpression of FTD-associated tau variants to model tau aggregation in hiPSC-derived neurons [46, 47]. The slow spontaneous aggregation rate of these variants, however, in combination with extracellular tau as a possible means of disease propagation [48, 49], led to the development of tau seeding models. In this paradigm, tau aggregation can be induced by exposure to pre-formed tau aggregates (seeds) in rodent [50] and human [46] cell models as well as in transgenic mice [51]. Tau seeding was also implemented for hiPSC-derived neurons embedded in a 7.5-fold lower concentration ECM than used in our model [47], which in our hands was not sufficient to generate a 3D structure. More importantly, seeds can induce cell non-autonomous responses independently of intraneuronal tau pathology [52, 53]. Since this may be a confounding factor for the investigation of how tau pathology within neurons affects their interaction with astrocytes in the 3D culture, we adopted a novel approach for seed-independent tau pathology using double mutant FTDtau1+2. This approach previously enabled seed-independent tau pathology in rodent brain slice cultures [29] and a 2D human neuron/astrocyte co-culture [8]. In agreement, neuron-specific expression of FTDtau¹⁺² induces intraneuronal MC1-positive tau pathology within 4 weeks in the 3D co-culture. In addition, the pathological inclusions label positive for the thiophene dye HS169 [30], demonstrating the accumulation of stacked beta-sheets that are typical for inclusions in the human brain. Altogether, these data demonstrate the utility of the 3D co-culture to model intraneuronal tau pathology that does not require seeds.

Our system is highly versatile and can easily be adapted to address more specific research questions. For example, the protocol contains a 24-h time window for manipulation with lentiviral transduction before co-culture. This allows the introduction of cell-specific manipulations like overexpression or knock-down of (e.g. disease-associated) proteins without the need for cell type-specific promoters. Moreover, the model can be used for the study of neuron-neuron interactions. As a proof-of-concept, we mixed two neuron populations overexpressing distinctly labelled FTDtau variants, demonstrating the specificity of FTDtau1+2 for the formation of tau pathology in the same 3D culture. In tau transgenic mice it was demonstrated that tau pathology can propagate to neurons as well as astrocytes that did not express the transgene [54]. Since in the mixed 3D co-culture of FTDtau¹- and FTDtau¹⁺²-expressing neurons we only observe tau pathology in the latter population, our data suggests that tau pathology in the population with FTDtau¹⁺² did not propagate within the timeline of our experiments.

Importantly, despite the close neuron-astrocyte contact, we do not observe transfer of tau (pathology) between neurons and astrocytes in the 3D culture. Moreover, we do not observe overt effects of tau pathology on neurons or astrocytes in the 3D co-culture. This makes the model suitable to study early events in tau pathogenesis that precede neurodegeneration.

A limitation of this study is the qualitative assessment of cellular phenotypes. More subtle and/or cell nonautonomous effects of tau pathology should be addressed in a quantitative manner in future studies. To this end, future efforts should prioritize the development of advanced cell segmentation tools for high-content imaging in 3D. Readouts can include the loss of synapses and neurites, the extent of tau pathology, the transfer of tau between neurons and astrocytes, as well as astrocytic morphology and (inflammatory) marker expression. In addition, single-cell omics techniques can be used to assess the cell-specific neuronal and astrocytic transcriptional responses to intraneuronal tau pathology. The feasibility of this approach to assess cell type-specific disease pathways in 3D cultures was recently demonstrated in a single-cell RNA sequencing study of iPSC-derived organoids, showing that tau mutations lead to cholesterol dysregulation specifically in astrocytes [55], but tau pathology was not studied. Since preliminary data indicate that our 3D co-cultures are also stable at 8 weeks, both these quantitative approaches may be complemented with prolonged culturing.

Our model does not include all cellular components: also other brain cells and in particular microglia play a crucial role in the pathogenesis and progression of tau pathology [56, 57]. Indeed, microglial activation coincides with tau pathology in patients [58-60] and rodent studies have demonstrated a variety of microglial processes that contribute to the progression of the disease, including the phagocytosis of neurons, synapses and tau, as well as the secretion of pro-inflammatory cytokines [61–66]. Rodent models have strongly contributed to our understanding of the (micro)glial component of tauopathies. Although gene expression patterns are remarkably similar between rodent and human microglia, some striking species-specific differences exist in microglial immune regulation [67] and aging [68]. In addition, the cellular and anatomical complexity of in vivo models may limit further mechanistic insight and cell type-specific intervention in the communication between neurons, astrocytes and microglia. On the other hand, 2D in vitro monocultures lack the in vivo complexity to accurately model the physiological and pathological responses of neural cells in general, and microglia in particular [67, 69] similar to astrocytes, as mentioned above. For ß-amyloid formation, a 3D neuron/astrocyte/microglia tri-culture has been developed using immortalized microglia [28], but no 3D tri-culture of tau aggregation exists yet. The 3D co-culture model of tau pathology we present here can be further developed to a 3D tri-culture by incorporating microglia. This approach would create the possibility to investigate interactions between microglia, astrocytes and neurons with tau pathology, and their role in neuroinflammation and neurodegeneration.

Conclusions

We have developed a 3D human co-culture model with synapse formation, neuronal activity and abundant neuron-astrocyte interactions. Using this model, we are the first to introduce intraneuronal tau pathology in a 3D co-culture in a cell-specific and seed-independent manner. This model therefore provides a suitable miniaturized platform to study cell (non-) autonomous disease mechanisms of tau pathogenesis in a physiologically relevant 3D environment.

Methods

Stem cell and primary human astrocyte cultures

A previously generated hiPSC line stably expressing a doxycyclin-inducible rtTA/Ngn2 [14] was maintained feeder-free on Vitronectin (Stem Cell Technologies) coated plates in TeSR-E8 medium (Stem Cell Technologies; 1X TeSR-E8 supplement (Stem Cell Technologies), 50 U/mL Pen/strep (Gibco)) supplemented with G418 (50 μg/mL; Sigma) and puromycin (0.5 μg/mL; Sigma) at 37 °C/5% CO₂ under hypoxic conditions. Colonies were fed daily and double-volume feeding allowed for weekend-free culturing. Once or twice a week, colonies were passaged by non-enzymatic dissociation with Gentle Cell Dissociation Reagent (Stem Cell Technologies) and replated in the presence of Rho-associated, coiled-coil containing protein kinase (ROCK) inhibitor (RI; 10 µM; Selleckchem) for 24 h to promote cell survival. To prevent genomic instability and senescence by prolonged passaging of hiPSCs, cryopreserved cells were thawed 2 weeks prior to use once cells reached passage number 25.

Primary human fetal astrocytes (ScienCell #1800) were cultured in Astrocyte Medium (ScienCell; 1X Growth Supplement, 1X FBS, 100 U/mL Pen/strep (all ScienCell)) on Geltrex (1:100 in DMEM/F12 (both Gibco) coated plates and half of the medium was refreshed twice a week. Astrocytes were passaged once a week by enzymatic dissociation with Accutase (Sigma). Cultures were kept at 37 °C/5% CO₂ at physiological oxygen conditions and used until passage number 7. To eliminate variability in astrocyte physiology due to genetic background of the donor, all experiments in this study were conducted with cryopreserved cells that were expanded from one vial.

Human co-culture and neuronal differentiation

The materials and standard operating procedures to generate the 3D human co-cultures are described in Supplementary Document 1. Briefly, two days prior to 3D co-culture (i.e. day -2), hiPSCs were non-enzymatically dissociated with Gentle Cell Dissociation Reagent and plated as single cells (80.000 cells/well) in 6-well culture plates (VWR) precoated with Geltrex (1:100 in DMEM/F12) in TesR-E8 medium with NT3 (10 ng/mL; Peprotech), BDNF (10 ng/mL; Peprotech), RI (10 µM) and doxycycline (2 µg/mL; Sigma) for induction of the Ngn2 transgene. After one day (i.e. day -1), medium was changed to DMEM/F12 (Gibco; 1X N2 supplement, 1X NEAA, 100 U/mL Pen/strep (all Gibco)) with NT3 (10 ng/mL), BDNF (10 ng/mL) and doxycycline (2 µg/ mL). At the day of co-culture (i.e. day 0), Ngn2 neural progenitors and primary human astrocytes were treated with Accutase for detachment. Prior to detachment, Accutase was aspirated and cells were dissociated separately in Neurobasal Medium (Gibco; 1X B27 + VitA, 1X Glutamax, 100 U/mL Pen/strep (all Gibco)) with NT3 (10 ng/mL), BDNF (10 ng/mL), RI (10 µM) and doxycyclin (2 µg/mL). After centrifugation the supernatant was removed and the co-culture cell suspension (30.000 neural precursors with 5.000 human astrocytes in 25 μL per well) was mixed with undiluted Geltrex (50% v/v final concentration) and kept on ice to prevent polymerization. Co-culture suspensions were plated in black 96-well culture plates (Greiner #655,090) and, once a week, half (week 1 and 2) or two-third (week 3) of the medium was refreshed until week 4. Only the inner 32 wells of the 96-well plate were used for culturing and empty wells were filled with sterile PBS to minimize evaporation of the culture medium. Throughout the neuronal differentiation process, co-cultures were kept at 37 °C/5% CO₂ at physiological oxygen conditions. See Supplementary Document 1 for a step-by-step culture procedure.

Lentiviral transduction

Lentiviral transduction of Ngn2-neural progenitors was performed one day prior to co-culture (i.e. day -1) with medium change to DMEM/F12 containing 4 μ L/mL viral solution. At the day of co-culture (i.e. day 0), 24 h after transduction, cells were washed once for 30 min with pre-warmed Neurobasal Medium containing all factors to remove lentiviral particles and thus prevent transduction of astrocytes during subsequent plating. In a 2nd generation lentiviral backbone, the expression of 2N4R FTDtau with an in-frame EGFP or mCherry C-terminal fluorescent tag was driven by a cytomegalovirus promoter. The following FTDtau mutants were used: P301L (FTDtau¹) and P301L+S320F (FTDtau¹⁺²). Lentiviral particles were generated as previously described [70].

Immunofluorescence staining

Culture media was removed and cells were washed with PBS and fixed with paraformaldehyde (PFA), consisting of 30 min incubation at room temperature (RT) in 4% PFA (Electron Microscopy Sciences) in PBS (pH 7.4). Fixed cultures were washed once with PBS for 30 min and immediately processed for immunostaining or stored at 4 °C. For immunostaining, cells were blocked and permeabilized for 3 h in blocking solution (5% normal goat serum (Thermo Fisher) + 0.1% bovine serum albumin and 0.3% Triton X-100 (both Sigma) in PBS (pH 7.4)), and subsequently incubated with primary antibodies diluted in blocking solution for 1 h at RT followed by overnight at 4 °C. After washing three times in PBS, cells were incubated overnight at 4 °C with Alexa Fluor-conjugated secondary antibodies (488, 546 or 568 and 647) (Invitrogen) diluted 1:1000 in blocking solution. HS169 (1:500, a kind gift of Peter Nilsson) was used during secondary antibody incubation to label stacked beta-sheets [30]. Cells were washed three times in PBS for 15 min at RT, and cell nuclei were stained with 4',6-diamidino-2-phenylindole (DAPI) (2.5 µg/mL, Sigma) in PBS of the second wash. Confocal images were acquired with a Nikon Eclipse Ti confocal microscope equipped with a $10 \times air$ and $40 \times oil$ immersion objective and controlled by NisElements 4.30 software (Nikon). Z-stacks with a step size between 1 μm (for stacks $< 15 \mu m$) and 5 μm (for stacks $> 15 \mu m$) were obtained and each parameter, including laser power and detector gain, was kept constant within each experiment for valid comparison. The following primary antibodies were used: CD44 (1:100, DSHB, H4C4), Cy3-conjugated NeuroChrom (1:500, Merck, ABN2300C3), GFAP (1:1000, DAKO, Z0334), MAP2 (1:2500, Millipore, AB5543), MC1 (1:500, kind gift of Peter Davies), NF200 (1:200, Sigma, N4142), SMI312 (1:1000, Eurogentec, 312P-050), PSD95 (1:200, Abcam, AB2723), Synaptophysin 1 (1:500, Synaptic Systems, 101,004) and Tuj1 (1:1000, R&D systems, MAB1195). Maximum intensity projected Z-stacks are shown. ImageJ software (National Institute of Health) was used to adapt images for publication. Immunofluorescence data are supported by at least 3 independent experiments.

Calcium imaging

To define spontaneous calcium transients, cells were incubated with 1 μ M Fluo-5F-AM (Molecular Probes, F14222) in complete culture medium for 30 min at 37 °C/5% CO₂. High-speed confocal imaging was performed in resonant mode using the Nikon Eclipse Ti confocal microscope, equipped with a 10 × air objective and controlled by NisElements 4.30 software (Nikon). Images (512 × 256 pixels at 2 × scanner zoom) were acquired for 10 min at 30 Hz on a 37 °C heated stage. Raw intensity

values were extracted in ImageJ (National Institute of Health) for each cell source, and the lowest and highest fluorescent signal was scaled between 0 and 1 respectively. Calcium data are supported by 3 independent experiments.

Electron microscopy

Culture media was removed and cells were washed once with PBS, followed by fixation with 2.5% glutaraldehyde (Merck) in 0.1 M cacodylate buffer, pH7.4, overnight at 4 °C. Samples were washed three times with 0.1 M cacodylate buffer and post-fixed by 1% OsO₄/1% Ru(CN)₆ for 90 min at RT. After washing with 0.1 M cacodylate buffer, a series of increasing ethanol concentrations (30%-100%) were used to dehydrate the samples. Samples were embedded in EPON and polymerized for 48 h at 65 °C. The EPON-embedded samples were removed from the substrate by breaking the 96-well plate in a vise. Regions with highest cell density were determined by light microscopy and ultrathin sections of (80 nm) were cut using an Ultracut UCT microtome (Leica). Sections were collected on single-slot, formvar-coated copper grids and counterstained with 0.5% uranyl acetate for 25 min and lead cytrate for 3 min at RT in an automated section stainer (Leica). Images were captured on a JEOL 1010 transmission EM at 60 kV with a side-mounted CCD camera (Modera, EMSIS) controlled by the iTEM software (EMSIS GmbH, Germany) at 40 K or 80 K magnification. EM data are supported by 2 independent experiments.

Abbreviations

2D Two-dimensional
 3D Three-dimensional
 AD Alzheimer's disease
 ECM Extracellular matrix
 EM Electron microscopy
 FTD Frontotemporal dementia

hiPSCs Human induced pluripotent stem cells

Ngn2 Neurogenin 2

PSD95 Post-synaptic density protein 95

SYP1 Synaptophysin 1

Supplementary Information

The online version contains supplementary material available at https://doi.org/10.1186/s12575-023-00194-2.

Additional file 1: Supplementary Figure 1. Annotation of calcium activity by immunofluorescence. A 4-week old 3D human neuron/ astrocyte co-culture was loaded with Fluo5-AM for calcium imaging and subsequently labeled with β-3-tubulin and GFAP to assess the neuronal/ astrocytic source of the calcium signal. (A) A 10-minute maximum fluorescence intensity projection of Fluo5-AM. The heatmap corresponds to the calcium concentration. The numbers correspond to the cells in B and calcium traces in C. See Supplementary Video 2 for calcium signals over time. (B) After calcium imaging, immunofluorescent staining was performed for neurons (β-3-tubulin, green) and astrocytes (GFAP, magenta), and nuclei were visualized using DAPI (greyscale). The numbers correspond to the

individual calcium signals in A and traces in C. (C) The graph shows the Fluo5-AM fluorescence intensity of the cells numbered in B cell over 10 minutes, with the lowest and highest value set to 0 and 1 respectively. Intracellular calcium transients were annotated using immunofluorescence for neurons (green) or astrocytes (purple). For two traces (cell 10 and 12), cellular sources were unidentified (black) due to the large overlap of β -3-tubulin and GFAP. Note that the calcium signal of one neuron (cell 9) remains high, possibly by cell death during imaging, as well as that the signal in one astrocyte (cell 15) slowly increases, possibly by passive uptake of the Fluo5-AM dye.

Additional file 2: Supplementary Document 1. Detailed list of reagents and standard operating procedures to generate the 3D human neuron/astrocyte co-cultures.

Additional file 3: Supplementary Video 1. Astrocytes enwrap neuronal cell bodies and neurites. Z-stack confocal video of neurons and astrocytes in a 112 µm-thick 3D co-culture displayed in Figure 2B. Immunostaining was performed for neurons (NeuroChrom, green) and astrocytes (GFAP, magenta).

Additional file 4: Supplementary Video 2. Calcium activity in 3D neuron/astrocyte co-culture. A 4-week old 3D human neuron/astrocyte co-culture was loaded with Fluo5-AM for calcium imaging. Confocal video showing Fluo5-AM fluorescence over time for 10 minutes in the 3D co-culture shown in Figure 3D and Supplementary Figure 1A.

Acknowledgements

We thank Nicole Breeuwsma and Lisa Gasparotto for advice and expert technical assistance on hiPSC- cultures; Rien Dekker for acquisition and help with interpretation of EM data; Jurjen Broeke for technical assistance on calcium imaging; Robbert Zalm, Ingrid Saarloos and Joost Hoetjes for preparation of constructs and viral particles. Also, we thank Peter Davies (Albert Einstein College of Medicine, USA) for providing the MC1 antibody, and Peter Nilsson (Linköping University, Sweden) for providing the HS169 dye. For critical reading of the manuscript and/or stimulating discussions we thank the members of the Molecular Neurodegeneration group (Dept. of Functional Genomics, Vrije Universiteit Amsterdam, Amsterdam) and the Stem Cell Biology group (Dept. of Complex Trait Genetics, Vrije Universiteit Amsterdam, Amsterdam).

Authors' contributions

K.L.B., L.S.P., V.M.H. and W.S. conceptualized the project. K.L.B., C.S. P.C-S. optimized the culture procedures and generated/processed the 3D cultures. K.L.B. conducted and analyzed the experiments. J.R.T. v W. acquired and interpreted EM data. K.L.B. and W.S. interpreted the data and wrote the manuscript with input from all authors. The author(s) read and approved the final manuscript.

Funding

This study was supported by ZonMW and Stichting Proefdiervrij (MKMD #114022506 to WS) and co-funded by the PPP Allowance made available by Health ~ Holland, Top Sector Life Sciences & Health, to stimulate public–private partnerships (#LSHM17014 to VMH and LSHM18024 to WS).

Availability of data and materials

The data that support the findings of this study are available from the corresponding author upon reasonable request.

Declarations

Ethics approval and consent to participate

Not applicable.

Consent for publication

Not applicable.

Competing interests

C.S. and L.S.P. are employees of Crown Bioscience B.V. (formerly OcellO B.V.), Leiden, The Netherlands.

Received: 2 September 2022 Accepted: 17 February 2023 Published online: 23 February 2023

References

- 1. Leyns CEG, Holtzman DM. Glial contributions to neurodegeneration in tauopathies. Mol Neurodegener. 2017;12:1–16.
- 2. Sidoryk-Węgrzynowicz M, Struzyńska L. Astroglial contribution to taudependent neurodegeneration. Biochem J. 2019;476:3493–504.
- Sidoryk-Wegrzynowicz M, Gerber YN, Ries M, Sastre M, Tolkovsky AM, Spillantini MG. Astrocytes in mouse models of tauopathies acquire early deficits and lose neurosupportive functions. Acta Neuropathol Commun. 2017:5:89.
- Jiwaji Z, Tiwari SS, Avilés-Reyes RX, Hooley M, Hampton D, Torvell M, et al. Reactive astrocytes acquire neuroprotective as well as deleterious signatures in response to Tau and Aß pathology. Nat Commun. 2022;13:1–23.
- Reid MJ, Beltran-Lobo P, Johnson L, Perez-Nievas BG, Noble W. Astrocytes in Tauopathies. Front Neurol. 2020;11:572850.
- Drummond E, Wisniewski T. Alzheimer's disease: experimental models and reality. Acta Neuropathol. 2017;133:155–75.
- Sidoryk-Wegrzynowicz M, Strużyńska L. Dysfunctional glia: contributors to neurodegenerative disorders. Neural Regen Res. 2021;2:218–22.
- Batenburg KL, Kasri NN, Heine VM, Scheper W. Intraneuronal tau aggregation induces the integrated stress response in astrocytes. J Mol Cell Biol. 2022. https://doi.org/10.1093/JMCB/MJAC071.
- Centeno EGZ, Cimarosti H, Bithell A. 2D versus 3D human induced pluripotent stem cell-derived cultures for neurodegenerative disease modelling. Mol Neurodegener. 2018;13(1):27.
- Watson PMD, Kavanagh E, Allenby G, Vassey M. Bioengineered 3D glial cell culture systems and applications for neurodegeneration and Neuroinflammation. SLAS Discov. 2017;22:583–601.
- Krencik R, Seo K, van Asperen JV, Basu N, Cvetkovic C, Barlas S, et al. Systematic three-dimensional COCULTURE rapidly recapitulates interactions between human neurons and astrocytes. Stem Cell Rep. 2017;9:1745–53.
- Pasca AM, Sloan SA, Clarke LE, Tian Y, Makinson CD, Huber N, et al. Functional cortical neurons and astrocytes from human pluripotent stem cells in 3D culture. Nat Methods. 2015;12:671–8.
- 13. Lu X, Yang J, Xiang Y. Modeling human neurodevelopmental diseases with brain organoids. Cell Regen. 2022;11:1–13.
- Frega M, Van Gestel SHC, Linda K, Van Der Raadt J, Keller J, Van Rhijn JR, et al. Rapid neuronal differentiation of induced pluripotent stem cells for measuring network activity on micro-electrode arrays. J Vis Exp. 2017;2017:1–10.
- Zhang Y, Pak CH, Han Y, Ahlenius H, Zhang Z, Chanda S, et al. Rapid single-step induction of functional neurons from human pluripotent stem cells. Neuron. 2013;78:785–98.
- Thoma EC, Wischmeyer E, Offen N, Maurus K, Sirén AL, Schartl M, et al. Ectopic expression of neurogenin 2 alone is sufficient to induce differentiation of embryonic stem cells into mature neurons. PLoS ONE. 2012;7:e38651.
- 17. Zhou B, Zuo YX, Jiang RT. Astrocyte morphology: Diversity, plasticity, and role in neurological diseases. CNS Neurosci Ther. 2019;25:665–73.
- Middeldorp J, Hol EM. GFAP in health and disease. Prog Neurobiol. 2011;93:421–43.
- Sosunov AA, Wu X, Tsankova NM, Guilfoyle E, McKhann GM, Goldman JE. Phenotypic heterogeneity and plasticity of isocortical and hippocampal astrocytes in the human brain. J Neurosci. 2014;34:2285–98.
- 20. Wevers NR, van Vught R, Wilschut KJ, Nicolas A, Chiang C, Lanz HL, et al. High-throughput compound evaluation on 3D networks of neurons and glia in a microfluidic platform. Sci Rep. 2016;6:1–10.
- Papadimitriou C, Celikkaya H, Cosacak MI, Mashkaryan V, Bray L, Bhattarai P, et al. 3D Culture Method for Alzheimer's Disease Modeling Reveals Interleukin-4 Rescues Aβ42-Induced Loss of Human Neural Stem Cell Plasticity. Dev Cell. 2018;46:85-101.e8.
- Oberheim NA, Takano T, Han X, He W, Lin JHC, Wang F, et al. Uniquely hominid features of adult human astrocytes. J Neurosci. 2009;29:3276–87.
- Allen NJ, Eroglu C. Cell Biology of Astrocyte-Synapse Interactions. Neuron. 2017;96:697–708.
- 24. Meijer M, Rehbach K, Brunner JW, Classen JA, Lammertse HCA, van Linge LA, et al. A single-cell model for synaptic transmission and

- plasticity in human iPSC-derived neurons. Cell Rep. 2019;27:2199-2211. e6.
- 25. Verstraelen P, Garcia-Diaz Barriga G, Verschuuren M, Asselbergh B, Nuydens R, Larsen PH, et al. Systematic Quantification of Synapses in Primary Neuronal Culture. Science. 2020;23:101542.
- Bélanger M, Allaman I, Magistretti PJ. Brain energy metabolism: focus on astrocyte-neuron metabolic cooperation. Cell Metab. 2011;14:724–38.
- 27. Grienberger C, Konnerth A. Imaging calcium in neurons. Neuron. 2012;73:862–85.
- Park J, Wetzel I, Marriott I, Dréau D, D'Avanzo C, Kim DY, et al. A 3D human triculture system modeling neurodegeneration and neuroinflammation in Alzheimer's disease. Nat Neurosci. 2018;21:941–51.
- Croft CL, Cruz PE, Ryu DH, Ceballos-Diaz C, Strang KH, Woody BM, et al. rAAV-based brain slice culture models of Alzheimer's and Parkinson's disease inclusion pathologies. J Exp Med. 2019;216:539–55.
- Shirani H, Linares M, Sigurdson CJ, Lindgren M, Norman P, Nilsson KPR.
 A palette of fluorescent thiophene-based ligands for the identification of protein aggregates. Chem Eur J. 2015;21:15133–7.
- 31. Andrews MG, Kriegstein AR. Challenges of Organoid Research. Annu Rev Neurosci. 2022;45:23–39.
- 32. de Leeuw SM, Davaz S, Wanner D, Milleret V, Ehrbar M, Gietl A, et al. Increased maturation of iPSC-derived neurons in a hydrogel-based 3D culture. J Neurosci Methods. 2021;360:109254.
- Balasubramanian S, Packard JA, Leach JB, Powell EM. Three-dimensional environment sustains morphological heterogeneity and promotes phenotypic progression during astrocyte development. Tissue Eng Part A. 2016;22:885–98.
- Placone AL, McGuiggan PM, Bergles DE, Guerrero-Cazares H, Quiñones-Hinojosa A, Searson PC. Human astrocytes develop physiological morphology and remain quiescent in a novel 3D matrix. Biomaterials. 2015;42:134–43.
- Stogsdill JA, Ramirez J, Liu D, Kim YH, Baldwin KT, Enustun E, et al. Astrocytic neuroligins control astrocyte morphogenesis and synaptogenesis. Nature. 2017;551:192–7.
- 36. Hasel P, Dando O, Jiwaji Z, Baxter P, Todd AC, Heron S, et al. Neurons and neuronal activity control gene expression in astrocytes to regulate their development and metabolism. Nat Commun. 2017;8 May 1:1–18.
- Pannasch U, Freche D, Dallérac G, Ghézali G, Escartin C, Ezan P, et al. Connexin 30 sets synaptic strength by controlling astroglial synapse invasion. Nat Neurosci. 2014;17:549–58.
- Verkhratsky A, Nedergaard M. Physiology of astroglia. Physiol Rev. 2018;98:239–389.
- Brennand K, Savas JN, Kim Y, Tran N, Simone A, Hashimoto-Torii K, et al. Phenotypic differences in hiPSC NPCs derived from patients with schizophrenia. Mol Psychiatry. 2015;20:361–8.
- Sposito T, Preza E, Mahoney CJ, Setó-Salvia N, Ryan NS, Morris HR, et al. Developmental regulation of tau splicing is disrupted in stem cellderived neurons from frontotemporal dementia patients with the 10 + 16 splice-site mutation in MAPT. Hum Mol Genet. 2015;24:5260–9.
- 41. Iovino M, Agathou S, González-Rueda A, Del Castillo V-H, Borroni B, Alberici A, et al. Early maturation and distinct tau pathology in induced pluripotent stem cell-derived neurons from patients with MAPT mutations. Brain. 2015;138:3345–59.
- Verheyen A, Diels A, Reumers J, Van Hoorde K, Van den Wyngaert I, van Outryve DC, et al. Genetically engineered iPSC-Derived FTDP-17 MAPT neurons display mutation-specific neurodegenerative and neurodevelopmental phenotypes. Stem Cell Rep. 2018;11:363–79.
- Paonessa F, Evans LD, Solanki R, Larrieu D, Wray S, Hardy J, et al. Microtubules deform the nuclear membrane and disrupt nucleocytoplasmic transport in tau-mediated frontotemporal dementia. Cell Rep. 2019;26:582-593.e5.
- 44. Ehrlich M, Hallmann AL, Reinhardt P, Araúzo-Bravo MJ, Korr S, Röpke A, et al. Distinct neurodegenerative changes in an induced pluripotent stem cell model of frontotemporal dementia linked to mutant TAU protein. Stem Cell Rep. 2015;5:83–96.
- Fong H, Wang C, Knoferle J, Walker D, Balestra ME, Tong LM, et al. Genetic correction of tauopathy phenotypes in neurons derived from human induced pluripotent stem cells. Stem Cell Rep. 2013;1:226–34.

- Verheyen A, Diels A, Dijkmans J, Oyelami T, Meneghello G, Mertens L, et al. Using human iPSC-derived neurons to model TAU aggregation. PLoS ONE. 2015;10:1–15.
- 47. Medda X, Mertens L, Versweyveld S, Diels A, Barnham L, Bretteville A, et al. Development of a scalable, high-throughput-compatible assay to detect tau aggregates using iPSC-derived cortical neurons maintained in a three-dimensional culture format. J Biomol Screen. 2016;21:804–15.
- Yamada K. Extracellular tau and its potential role in the propagation of tau pathology. Front Neurosci. 2017;11 Nov:667.
- 49. Frost B, Jacks RL, Diamond Ml. Propagation of Tau misfolding from the outside to the inside of a cell. J Biol Chem. 2009;284:12845–52.
- Wiersma VI, van Ziel AM, Vazquez-Sanchez S, Nölle A, Berenjeno-Correa E, Bonaterra-Pastra A, et al. Granulovacuolar degeneration bodies are neuron-selective lysosomal structures induced by intracellular tau pathology. Acta Neuropathol. 2019;138:943

 –70.
- 51. Peeraer E, Bottelbergs A, van Kolen K, Stancu IC, Vasconcelos B, Mahieu M, et al. Intracerebral injection of preformed synthetic tau fibrils initiates widespread tauopathy and neuronal loss in the brains of tau transgenic mice. Neurobiol Dis. 2015;73:83–95.
- Wang P, Ye Y. Filamentous recombinant human Tau activates primary astrocytes via an integrin receptor complex. Nat Commun. 2021;12:1–14.
- Ungerleider K, Beck J, Lissa D, Turnquist C, Horikawa I, Harris BT, et al. Astrocyte senescence and SASP in neurodegeneration: tau joins the loop. Cell Cycle. 2021;20:752–64.
- de Calignon A, Polydoro M, Suárez-Calvet M, William C, Adamowicz DH, Kopeikina KJ, et al. Propagation of tau pathology in a model of early Alzheimer's disease. Neuron. 2012;73:685–97.
- Glasauer SMK, Goderie SK, Rauch JN, Guzman E, Audouard M, Bertucci T, et al. Human tau mutations in cerebral organoids induce a progressive dyshomeostasis of cholesterol. Stem Cell Rep. 2022;17:2127–40.
- 56. Didonna A. Tau at the interface between neurodegeneration and neuro-inflammation. Genes Immun. 2020;21:288–300.
- Salter MW, Stevens B. Microglia emerge as central players in brain disease. Nat Med. 2017;23:1018–27.
- Sheffield LG, Marquis JG, Berman NEJ. Regional distribution of cortical microglia parallels that of neurofibrillary tangles in Alzheimer's disease. Neurosci Lett. 2000;285:165–8.
- Serrano-Pozo A, Mielke ML, Gómez-Isla T, Betensky RA, Growdon JH, Frosch MP, et al. Reactive glia not only associates with plaques but also parallels tangles in Alzheimer's disease. Am J Pathol. 2011;179:1373–84.
- Pascoal TA, Benedet AL, Ashton NJ, Kang MS, Therriault J, Chamoun M, et al. Microglial activation and tau propagate jointly across Braak stages. Nat Med. 2021;2021(27):1592–9.
- Asai H, Ikezu S, Tsunoda S, Medalla M, Luebke J, Haydar T, et al. Depletion of microglia and inhibition of exosome synthesis halt tau propagation. Nat Neurosci. 2015;18:1584–93.
- 62. Brelstaff J, Tolkovsky AM, Ghetti B, Goedert M, Spillantini MG. Living neurons with tau filaments aberrantly expose phosphatidylserine and are phagocytosed by microglia. Cell Rep. 2018;24:1939-1948.e4.
- Brelstaff JH, Mason M, Katsinelos T, McEwan WA, Ghetti B, Tolkovsky AM, et al. Microglia become hypofunctional and release metalloproteases and tau seeds when phagocytosing live neurons with P301S tau aggregates. Sci Adv. 2021:7:eabq4980.
- Dejanovic B, Huntley MA, De Mazière A, Meilandt WJ, Wu T, Srinivasan K, et al. Changes in the synaptic proteome in Tauopathy and rescue of tauinduced synapse loss by C1q antibodies. Neuron. 2018;100:1322-1336.e7.
- Yoshiyama Y, Higuchi M, Zhang B, Huang SM, Iwata N, Saido TCC, et al. Synapse loss and microglial activation precede tangles in a P301S tauopathy mouse model. Neuron. 2007;53:337–51.
- Wang C, Fan L, Khawaja RR, Liu B, Zhan L, Kodama L, et al. Microglial NF-кB drives tau spreading and toxicity in a mouse model of tauopathy. Nat Commun. 2022;13:1–19.
- Gosselin D, Skola D, Coufal NG, Holtman IR, Schlachetzki JCM, Sajti E, et al. An environment-dependent transcriptional network specifies human microalia identity. Science. 2017;356:1248–59.
- Galatro TF, Holtman IR, Lerario AM, Vainchtein ID, Brouwer N, Sola PR, et al. Transcriptomic analysis of purified human cortical microglia reveals ageassociated changes. Nat Neurosci. 2017;20:1162–71.
- McFarland KN, Ceballos C, Rosario A, Ladd T, Moore B, Golde G, et al. Microglia show differential transcriptomic response to Aβ peptide aggregates ex vivo and in vivo. Life Sci Alliance. 2021;4:e202101108.

70. Naldini L, Blömer U, Gallay P, Ory D, Mulligan R, Gage FH, et al. In vivo gene delivery and stable transduction of nondividing cells by a lentiviral vector. Science. 1979;1996(272):263–7.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- $\bullet\,$ thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

